

## How to export your data from Zero-Click and import it into ALLY Home

1. Log in to your Zero-Click software and check the number of your **Profile ID** is identical to the **NHS ID** number you have entered in the **Personal Data** section of your ALLY Home account. If it is not identical you will need to complete step 2 and amend your **Profile ID** in Zero-Click to your **NHS ID**. Please note that data will only be imported into ALLY Home if both numbers are identical.
2. To amend your **Profile ID**, select **Profile > Setup Personal Info**, type in your **NHS ID** number in the **Profile ID** field and select **Save This Page**.
3. To begin the data export from Zero-Click select **File > File Maintenance**, the **File Maintenance Wizard** will then appear.
4. Within the **File Maintenance Wizard** select **Next > Export a profile > Next**. Select the profile you wish to export and select **Next**.
5. Select **Open format text file > Next**
6. The **File Maintenance Wizard** will ask you to select a file location. The file location field will already be prepopulated, and you will just need to select **Next**.
7. Select **Finish** and the data will automatically export into a csv file. The file will automatically save in the **Documents** folder on your PC.
8. Log in to ALLY Home and select **Personal Data > Import Zero-Click**. You will then be directed to your **Documents** folder.
9. Select the same csv file that is saved within your **Documents** folder and select **Open**.
10. ALLY Home will display a pop-up-message to say that the data transfer is complete. You will now be able to view your Zero-Click data in the **Data Management** section of ALLY Home.